

# mRNA Cap 2'-O-Methyltransferase




## M0366S



**2,000 units    50,000 U/ml    Lot: 0031403**  
**RECOMBINANT    Store at -20°C    Exp: 3/16**

**Description:** mRNA Cap 2'-O-Methyltransferase adds a methyl group at the 2'-O position of the first nucleotide adjacent to the cap structure at the 5' end of the RNA. The enzyme utilizes S-adenosylmethionine (SAM) as a methyl donor to methylate capped RNA (cap-0) resulting in a cap-1 structure.



The cap-1 structure has been reported to enhance mRNA translation efficiency (1) and hence may help improve expression in mRNA transfection and microinjection experiments.

mRNA Cap 2'-O-Methyltransferase specifically requires RNA with an m<sup>7</sup>GpppN cap as substrate. It cannot utilize RNA with pN, ppN, pppN or GpppN at the 5' end (2,3). Capped RNA may be prepared via *in vitro* transcription using cap analog or through enzymatic capping using the Vaccinia Capping Enzyme (NEB #M2080).

The reagents provided in this pack can be used to methylate up to 400 µg of capped RNA.

**Source:** An *E. coli* strain that carries the gene for the Vaccinia mRNA Cap 2'-O-Methyltransferase.

### Applications:

- 2'-O-methylation of capped mRNA for improved expression during microinjection and transfection experiments.

Supplied in: 20 mM Tris-HCl, pH 8.0, 100 mM NaCl, 1 mM DTT, 0.1 mM EDTA, 0.1% Triton X-100 and 50% glycerol.

### Reagents Supplied with Enzyme:

10X Capping Buffer    100 µl  
 SAM (32 mM)    100 µl

**Reaction Conditions:** 1X Capping Buffer.  
 Incubate at 37°C.

### 1X Capping Buffer:

50 mM Tris-HCl, pH 8.0  
 5 mM KCl  
 1 mM MgCl<sub>2</sub>  
 1 mM DTT

**Unit Definition:** One unit is defined as the amount of enzyme required to methylate 10 pmoles of 80 nt long capped RNA transcript in 1 hour at 37°C.

### Quality Control Assays

**RNase Assay:** Incubation of a 10 µl reaction containing 50 units of mRNA Cap 2'-O-Methyltransferase with 40 ng of 300-mer RNA transcript for 2 hours at 37°C resulted in less than 10% degradation of RNA as determined by denaturing PAGE analysis.

**Exonuclease Activity:** Incubation of a 50 µl reaction containing 50 units of mRNA Cap 2'-O-Methyltransferase with 1 µg of a mixture of single and double-stranded [<sup>3</sup>H] *E. coli* DNA for 4 hours at 37°C released < 0.5% of the total radioactivity.

**Endonuclease Activity:** Incubation of a 50 µl reaction containing 50 units of mRNA Cap 2'-O-Methyltransferase with 1 µg of φX174 RF I DNA for 4 hours at 37°C produced less than 10% conversion to RF II as determined by agarose gel electrophoresis.

**Protein Purity (SDS-PAGE):** mRNA Cap 2'-O-Methyltransferase is > 99% pure as determined by SDS PAGE analysis using Coomassie staining.

### Notes on Use

#### (READ PRIOR TO SETTING UP REACTION)

- RNA prepared using *in vitro* transcription and cap analog should be purified prior to use and resuspended in nuclease-free water. EDTA and salts should not be present in the solution.  
 mRNA Cap 2'-O-Methyltransferase may be directly added to a Vaccinia Capping System (NEB #M2080) reaction. RNA purification is not required in this case.
- Heating the RNA at 65°C for 5 minutes prior to incubation with the enzyme removes secondary structure on the 5' end of the transcript. Extend time to 10 minutes for transcripts with known highly structured 5' ends.
- SAM is unstable at pH 7-8, 37°C and should be mixed fresh prior to starting the reaction. We recommend determining how many reactions will be performed and diluting an aliquot of the 32 mM stock to 4 mM immediately before setting up the reactions. This "working stock" should be kept on ice to prevent degradation of SAM.

### Protocols

#### 2'-O-Methylation of Capped RNA

This protocol is designed to methylate up to 10 µg of capped RNA in a 20 µl reaction. Reaction size can be scaled up as needed.

- Combine capped RNA and nuclease-free water in a final volume of 16 µl. (Refer to step 1 in the notes on use).
- Heat at 65°C for 5 minutes (Refer to step 2 in the notes on use).
- Place tube on ice for 5 minutes.
- Add the following components in the order specified:
 

Denatured capped RNA (from above)	16.0 µl
10X Capping Buffer	2.0 µl
SAM (4 mM, dilute 32 mM stock to 4 mM)	1.0 µl
mRNA Cap 2'-O-Methyltransferase (50 U/µl)	1.0 µl
	20.0 µl

Note: Use of RNase Inhibitor is recommended to enhance stability of RNA in the reaction. Add 0.5 µl of RNase Inhibitor (e.g., Murine RNase Inhibitor NEB #M0314) during reaction set up. Subtract the additional volume from the amount of H<sub>2</sub>O used in the reaction.

- Incubate at 37°C for 60 minutes (For RNA less than 200 nt long increase incubation time to 2 hours)
- Proceed with purification of the RNA (if required) for downstream applications.

### One-Step Capping and 2'-O-Methylation

This protocol is designed to complete both capping and 2'-O-methylation in a single step. It involves incubating uncapped RNA with the Vaccinia Capping Enzyme (NEB #M2080, not included) and mRNA Cap 2'-O-Methyltransferase in the presence of GTP and SAM. The Vaccinia Capping Enzyme adds the cap at the 5' end of the RNA followed by 2'-O-methylation by the methyltransferase. This protocol can synthesize up to 10 µg of cap-1 RNA in a 20 µl reaction. Reaction size can be scaled up as needed.

- Combine uncapped RNA and nuclease-free water in a final volume of 14.0 µl. (Refer to step 1 in the notes on use).
- Heat at 65°C for 5 minutes (Refer to step 2 in the notes on use).
- Place tube on ice for 5 minutes
- Add the following components in the order specified:
 

Denatured RNA (from above)	14.0 µl
10X Capping Buffer	2.0 µl
GTP (10 mM)	1.0 µl
SAM (4 mM, dilute 32 mM stock to 4 mM)	1.0 µl
Vaccinia Capping Enzyme (10 U/µl)	1.0 µl
mRNA Cap 2'-O-Methyltransferase (50 U/µl)	1.0 µl
	20.0 µl

Note: Use of RNase Inhibitor is recommended to enhance stability of RNA in the reaction. Add 0.5 µl of RNase Inhibitor (e.g., Murine RNase Inhibitor NEB #M0314) during reaction set up. Subtract the additional volume from the amount of H<sub>2</sub>O used in the reaction.

- Incubate at 37°C for 60 minutes (For RNA less than 200 nt long increase incubation time to 2 hours)
- Proceed with purification of the RNA (if required) for downstream applications.

(see other side)

**References:**

1. Kuge, H. et al. (1998) *Nucleic Acids Res.* 26, 3208-3214.
2. Barbosa, E. and Moss, B. (1978) *J. Biol. Chem.* 253, 7698-7702.
3. Lockless, S.W. et al. (1998) *Biochemistry.* 37, 8564-8574.

**Companion Products**

## Vaccinia Capping System

M2080S 400 units

## RNase Inhibitor, Murine

M0314S 3,000 units

M0314L 15,000 units

## RNase Inhibitor, Human Placenta

M0307S 2,000 units

M0307L 10,000 units

## T7 High Yield RNA Synthesis Kit

E2040S 50 reactions

*E. coli* Poly(A) Polymerase

M0276S 100 units

M0276L 500 units

## Ribonucleotide Solution Set

N0450S 10  $\mu$ mol of eachN0450L 50  $\mu$ mol of each