
New England Biolabs Product Specification

<i>Product Name:</i>	<i>NEBNext[®] Multiplex Small RNA Library Prep Set for Illumina[®] (Set 2)</i>
<i>Catalog #:</i>	<i>E7580S/L</i>
<i>Kit Components:</i>	<i>NEBNext[®] 3' Ligation Reaction Buffer (E7301)</i> <i>NEBNext[®] 3' Ligation Enzyme Mix (E7288)</i> <i>NEBNext[®] 3' SR Adaptor for Illumina[®] (E7332)</i> <i>NEBNext[®] 5' Ligation Reaction Buffer (E7304)</i> <i>NEBNext[®] 5' Ligation Enzyme Mix (E7305)</i> <i>NEBNext[®] 5' SR Adaptor for Illumina[®] (E7328)</i> <i>NEBNext[®] SR RT Primer for Illumina[®] (E7333)</i> <i>NEBNext[®] First Strand Synthesis Reaction Buffer (E7334)</i> <i>ProtoScript[®] II Reverse Transcriptase (E7355)</i> <i>Murine RNase Inhibitor (E7308)</i> <i>LongAmp[®] Taq 2X Master Mix (E7309)</i> <i>NEBNext[®] SR Primer for Illumina[®] (E7310)</i> <i>Gel Loading Dye, Blue, 6X (E6138)</i> <i>Quick-Load[®] pBR322 DNA-MspI Digest (E7323)</i> <i>DNA Gel Elution Buffer (E7324)</i> <i>Linear Acrylamide (E7325)</i> <i>TE Buffer (E7326)</i> <i>Nuclease Free Water (E7327)</i> <i>NEBNext[®] Index 13 Primer for Illumina[®] (E7581)</i> <i>NEBNext[®] Index 14 Primer for Illumina[®] (E7582)</i> <i>NEBNext[®] Index 15 Primer for Illumina[®] (E7583)</i> <i>NEBNext[®] Index 16 Primer for Illumina[®] (E7584)</i> <i>NEBNext[®] Index 17 Primer for Illumina[®] (E7585)</i> <i>NEBNext[®] Index 18 Primer for Illumina[®] (E7586)</i> <i>NEBNext[®] Index 19 Primer for Illumina[®] (E7587)</i> <i>NEBNext[®] Index 20 Primer for Illumina[®] (E7588)</i> <i>NEBNext[®] Index 21 Primer for Illumina[®] (E7589)</i> <i>NEBNext[®] Index 22 Primer for Illumina[®] (E7590)</i> <i>NEBNext[®] Index 23 Primer for Illumina[®] (E7591)</i> <i>NEBNext[®] Index 24 Primer for Illumina[®] (E7592)</i>



New England Biolabs Product Specification

Shelf Life: 18 months
Storage Temp: -20°C
Specification Version: PS-E7580S/L v1.0
Effective Date: 15 Mar 2019

Assay Name/Specification (minimum release criteria)

<p>Functional Testing (Library Construction, Small RNA) - Each of the components is functionally validated and compared to the previous lot through construction of libraries made from commercially available human brain RNA using the kit's minimum and maximum input requirements. Libraries made from previous and current lots are sequenced on the same Illumina® flow cell and compared across various metrics including library yield and number of miRNAs identified.</p>
--

<p>* Individual Product Component Note - Standard Quality Control Tests are performed for each component included in NEBNext® Multiplex Small RNA Library Prep Set for Illumina® (Set 2) and meet the designated specifications.</p>



Date 15 Mar 2019

Derek Robinson
Director of Quality Control

