## **T4 DNA Ligase Reaction Buffer**







# **B0202S**

6.0 ml T4 DNA Ligase Reaction Buffer (10X) Store at -20°C Lot: 0021609 Exp: 9/19

**Description:** New England Biolabs supplies a 10X reaction buffer with all of its enzymes. At a 1X concentration this reaction buffer assures optimal activity of the enzyme.

#### 1X T4 DNA Ligase Reaction Buffer:

50 mM Tris-HCI 10 mM MgCl<sub>o</sub> 10 mM DTT 1 mM ATP pH 7.5 @ 25°C

Supplied with enzyme as a 10X concentrated stock

### **Quality Control Assays**

16-Hour Incubation: A 50 ul reaction containing this reaction buffer at a 1X concentration and 1 µg of HaelII digested  $\phi$ X174 RF I DNA incubated for 16 hours resulted in no detectable non-specific nuclease degradation.

Endonuclease Activity: Incubation of this reaction buffer at a 1X concentration with 1 μg φX174 RF I DNA for 4 hours at 37°C in 50 ul reactions resulted in less than 5% conversion to RF II.







NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans

CERTIFICATE OF ANALYSIS

## **T4 DNA Ligase Reaction Buffer**



1-800-632-7799 info@neb.com www.neb.com



# **B0202S**

6.0 ml T4 DNA Ligase Reaction Buffer (10X) Store at -20°C Lot: 0021609 Exp: 9/19

**Description:** New England Biolabs supplies a 10X reaction buffer with all of its enzymes. At a 1X concentration this reaction buffer assures optimal activity of the enzyme.

#### 1X T4 DNA Ligase Reaction Buffer:

50 mM Tris-HCI 10 mM MgCl<sub>o</sub> 10 mM DTT 1 mM ATP pH 7.5 @ 25°C Supplied with enzyme as a 10X concentrated stock

### **Quality Control Assays**

16-Hour Incubation: A 50 µl reaction containing this reaction buffer at a 1X concentration and 1 µg of HaelII digested \$\phi X174 RF I DNA incubated for 16 hours resulted in no detectable non-specific nuclease degradation.

Endonuclease Activity: Incubation of this reaction buffer at a 1X concentration with 1 μg φX174 RF I DNA for 4 hours at 37°C in 50 ul reactions resulted in less than 5% conversion to RF II.







NEW ENGLAND BIOLABS® is a registered trademark of New England

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans