

**β -Nicotinamide
adenine dinucleotide
(NAD⁺)**



1-800-632-7799
info@neb.com
www.neb.com



B9007S 008121014101

B9007S

0.2 ml 50 mM solution Lot: 0081210

Store at -20°C Exp: 10/14

Description: β -Nicotinamide adenine dinucleotide (NAD⁺) is included with the PreCR™ Repair Mix available from New England Biolabs. NAD⁺ is supplied as an aqueous 50 mM solution.

Use and Storage: For use in the PreCR Repair Mix, NAD⁺ should be diluted 1:100 to a final concentration of 500 μ M. NAD⁺ solution should be stored at -20°C, and the user should avoid repeated freeze/thaw cycles. For storage longer than 30 days, store at -80°C.

CERTIFICATE OF ANALYSIS

**β -Nicotinamide
adenine dinucleotide
(NAD⁺)**



1-800-632-7799
info@neb.com
www.neb.com



B9007S 008121014101

B9007S

0.2 ml 50 mM solution Lot: 0081210

Store at -20°C Exp: 10/14

Description: β -Nicotinamide adenine dinucleotide (NAD⁺) is included with the PreCR™ Repair Mix available from New England Biolabs. NAD⁺ is supplied as an aqueous 50 mM solution.

Use and Storage: For use in the PreCR Repair Mix, NAD⁺ should be diluted 1:100 to a final concentration of 500 μ M. NAD⁺ solution should be stored at -20°C, and the user should avoid repeated freeze/thaw cycles. For storage longer than 30 days, store at -80°C.

CERTIFICATE OF ANALYSIS