

SOC Outgrowth Medium



1-800-632-7799
info@neb.com
www.neb.com



B9020S 320140716071

B9020S

4 x 25 ml medium

Lot: 3201407

Store at 4°C

Exp: 7/16

Description: Add SOC to freshly transformed *E.coli* cells and incubate for up to 1 hour before plating on selective media.

1X SOC Outgrowth Medium:

2% Vegetable Peptone

0.5% Yeast Extract

10 mM NaCl

2.5 mM KCl

10 mM MgCl₂

10 mM MgSO₄

20 mM Glucose

Quality Assurance: SOC Medium is tested for sterility and optimal outgrowth of freshly transformed chemically competent and electrocompetent *E.coli*.


Storage Note: SOC Medium is shipped at 0°C or -80°C. Upon receipt, store at room temperature or 4°C.

CERTIFICATE OF ANALYSIS

SOC Outgrowth Medium



1-800-632-7799
info@neb.com
www.neb.com



B9020S 320140716071

B9020S

4 x 25 ml medium

Lot: 3201407

Store at 4°C

Exp: 7/16

Description: Add SOC to freshly transformed *E.coli* cells and incubate for up to 1 hour before plating on selective media.

1X SOC Outgrowth Medium:

2% Vegetable Peptone

0.5% Yeast Extract

10 mM NaCl

2.5 mM KCl

10 mM MgCl₂

10 mM MgSO₄

20 mM Glucose

Quality Assurance: SOC Medium is tested for sterility and optimal outgrowth of freshly transformed chemically competent and electrocompetent *E.coli*.

Storage Note: SOC Medium is shipped at 0°C or -80°C. Upon receipt, store at room temperature or 4°C.

CERTIFICATE OF ANALYSIS