HaeIII Methyltransferase



1-800-632-7799 info@neb.com www.neb.com







500 units 10.000 U/ml Lot: 0061403 RECOMBINANT Store at -20°C Exp: 3/15

Methylation Site:

5'... GGCC...3' 3'... CCGG...5' CH₂

Description: HaellI Methyltransferase modifies the internal cytosine residue (C5) in the sequence above.

Source: An *E. coli* strain that carries the cloned HaellI modification gene from *Haemophilus* aegyptius (ATCC 11116)

Supplied in: 50 mM KCl. 50 mM Tris-HCl (pH 7.5), 10 mM EDTA, 1 mM dithiothreitol. 200 µg/ml BSA and 50% glycerol.

Reagents Supplied with Enzyme:

10X Haelll Methyltransferase Reaction Buffer. 400X S-adenosylmethionine (32 mM).

Reaction Conditions: 1X Haelll

Methyltransferase Reaction Buffer, supplement

S-adenosylmethionine (supplied). Incubate at

1X Haelll Methyltransferase Reaction Buffer:

50 mM NaCl 50 mM Tris-HCL 10 mM dithiothreitol pH 8.5 @ 25°C

Protection Assav Conditions:

HaellI Methyltransferase is incubated with 1 μq of λ DNA in 10 μl of 1X Haelll Methyltransferase Reaction Buffer, supplemented

with 80 µM S-adenosylmethionine, for one hour at 37°C followed by 15 minutes at 65°C. The extent of protection is determined by the addition of 40 µl NEBuffer 2 and 10 units of HaellI restriction endonuclease. Incubation for 1 hour at 37°C is followed by analysis on an agarose gel.

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μ g λ DNA in 1 hour at 37°C in a total reaction volume of 10 ul against cleavage by HaelII restriction endonuclease.

Quality Control Assays

16-Hour Incubation: Incubation of 100 units with 1 μa of HindIII-diaested λ DNA in 50 μ I 1X NEBuffer 2 for 16 hours at 37°C resulted in no detectable endonuclease contamination.

Exonuclease Activity: Incubation of 100 units of HaeIII Methyltransferase with 1 µg sonicated ³H-DNA (10⁵ cpm/μg) for 4 hours at 37°C in 50 μl NEBuffer 2 [50 mM NaCl, 10 mM Tris-HCl (pH 7.9 @ 25°C), 10 mM MgCl₂, 1 mM dithiothreitol] released < 0.05% of the total radioactivity.

Endonuclease Contamination: Incubation of 100 units with 1 μg φX174 RF I DNA (4 hours, 37°C in 50 µl of NEBuffer 2) gave < 5% conversion to RF II.

Storage of SAM: S-adenosylmethionine (SAM) is stored at -20°C as a 32 mM solution dissolved in 0.005 M sulfuric acid and 10% ethanol. Under these conditions SAM is stable for up to 6 months. SAM is unstable at (pH 7.5), 37°C (1) and should be replenished in reactions incubated longer than 4 hours.

Methylation can be optimized by using fresh SAM.

Note: HaellI Methyltransferase protects DNA against cleavage by NotI.

Reference:

1. Hoffman, J.L. (1986) Biochemistry 25, 4444-4449.

Companion Product:

S-adenosylmethionine (SAM) #B9003S

0.5 ml

CERTIFICATE OF ANALYSIS

HaeIII Methyltransferase



1-800-632-7799 info@neb.com www.neb.com







Lot: 0061403

500 units 10.000 U/ml RECOMBINANT Store at -20°C Exp: 3/15

Methylation Site:

5'... GGCC...3' 3'... CCGG...5' CH,

Description: HaellI Methyltransferase modifies the internal cytosine residue (C5) in the sequence above.

Source: An *E. coli* strain that carries the cloned HaelII modification gene from *Haemophilus* aegyptius (ATCC 11116)

Supplied in: 50 mM KCl. 50 mM Tris-HCl (pH 7.5), 10 mM EDTA, 1 mM dithiothreitol. 200 µg/ml BSA and 50% glycerol.

Reagents Supplied with Enzyme:

10X Haelll Methyltransferase Reaction Buffer. 400X S-adenosylmethionine (32 mM).

Reaction Conditions: 1X Haelll

Methyltransferase Reaction Buffer, supplement with 80 uM S-adenosylmethionine (supplied). Incubate at

1X Haelli Methyltransferase Reaction Buffer:

50 mM NaCl 50 mM Tris-HCI 10 mM dithiothreitol pH 8.5 @ 25°C

Protection Assav Conditions:

HaellI Methyltransferase is incubated with 1 μq of λ DNA in 10 μl of 1X Haelll Methyltransferase Reaction Buffer, supplemented with 80 µM S-adenosylmethionine, for one hour at 37°C followed by 15 minutes at 65°C. The extent of protection is determined by the addition of 40 µl NEBuffer 2 and 10 units of HaelII restriction endonuclease. Incubation for 1 hour at 37°C is followed by analysis on an agarose gel.

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μ g λ DNA in 1 hour at 37°C in a total reaction volume of 10 µl against cleavage by HaelII restriction endonuclease.

Quality Control Assays

16-Hour Incubation: Incubation of 100 units with 1 μ g of HindIII-digested λ DNA in 50 μ I 1X NEBuffer 2 for 16 hours at 37°C resulted in no detectable endonuclease contamination.

Exonuclease Activity: Incubation of 100 units of HaellI Methyltransferase with 1 µg sonicated ³H-DNA (10⁵ cpm/μg) for 4 hours at 37°C in 50 μl NEBuffer 2 [50 mM NaCl, 10 mM Tris-HCl (pH 7.9 @ 25°C), 10 mM MgCl₂, 1 mM dithiothreitol] released < 0.05% of the total radioactivity.

Endonuclease Contamination: Incubation of 100 units with 1 ug ϕ X174 RF I DNA (4 hours. 37°C in 50 µl of NEBuffer 2) gave < 5% conversion to RF II.

Storage of SAM: S-adenosylmethionine (SAM) is stored at -20°C as a 32 mM solution dissolved in 0.005 M sulfuric acid and 10% ethanol. Under these conditions SAM is stable for up to 6 months. SAM is unstable at (pH 7.5), 37°C (1) and should be replenished in reactions incubated longer than 4 hours.

Methylation can be optimized by using fresh SAM.

Note: Haelli Methyltransferase protects DNA against cleavage by NotI.

Reference:

1. Hoffman, J.L. (1986) Biochemistry 25, 4444-4449.

Companion Product:

S-adenosylmethionine (SAM) #B9003S

0.5 ml