

## Antarctic Phosphatase



1-800-632-7799  
info@neb.com  
www.neb.com



M0289S 020131115111

# M0289S



**1,000 units**    **5,000 U/ml**    **Lot: 0201311**  
**RECOMBINANT**    **Store at -20°C**    **Exp: 11/15**

**Description:** Antarctic Phosphatase catalyzes the removal of 5' phosphate groups from DNA and RNA. Since phosphatase-treated fragments lack the 5' phosphoryl termini required by ligases, they cannot self-ligate (1). This property can be used to decrease the vector background in cloning strategies.

Heat Inactivated In 5 Minutes at 70°C

## Antarctic Phosphatase



1-800-632-7799  
info@neb.com  
www.neb.com



M0289S 020131115111

# M0289S



**1,000 units**    **5,000 U/ml**    **Lot: 0201311**  
**RECOMBINANT**    **Store at -20°C**    **Exp: 11/15**

**Description:** Antarctic Phosphatase catalyzes the removal of 5' phosphate groups from DNA and RNA. Since phosphatase-treated fragments lack the 5' phosphoryl termini required by ligases, they cannot self-ligate (1). This property can be used to decrease the vector background in cloning strategies.

Heat Inactivated In 5 Minutes at 70°C

**Source:** An *E. coli* strain that carries the TAB5 AP gene, originally cloned in plasmid pNI (2), recloned in plasmid pEGTAB7-4.1(3).

### Applications:

- Removing 5' phosphates from DNA, RNA, rNTPs and dNTPs
- Preparation of templates for 5' end labeling
- Prevention of recircularization of cloning vectors
- Dephosphorylation of proteins
- Removal of dNTPs and pyrophosphate from PCR reactions

Supplied in: 10 mM Tris-HCl (pH 7.4), 1 mM MgCl<sub>2</sub>, 0.01 mM ZnCl<sub>2</sub>, 1 mM DTT and 50% glycerol.

### Reagents Supplied with Enzyme:

10X Antarctic Phosphatase Reaction Buffer.

**Reaction Conditions:** 1X Antarctic Phosphatase Reaction Buffer. Incubate at 37°C.

### 1X Antarctic Phosphatase Reaction Buffer:

50 mM Bis-Tris-Propane HCl

1 mM MgCl<sub>2</sub>

0.1 mM ZnCl<sub>2</sub>

pH 6.0 @ 25°C

**Source:** An *E. coli* strain that carries the TAB5 AP gene, originally cloned in plasmid pNI (2), recloned in plasmid pEGTAB7-4.1(3).

### Applications:

- Removing 5' phosphates from DNA, RNA, rNTPs and dNTPs
- Preparation of templates for 5' end labeling
- Prevention of recircularization of cloning vectors
- Dephosphorylation of proteins
- Removal of dNTPs and pyrophosphate from PCR reactions

Supplied in: 10 mM Tris-HCl (pH 7.4), 1 mM MgCl<sub>2</sub>, 0.01 mM ZnCl<sub>2</sub>, 1 mM DTT and 50% glycerol.

### Reagents Supplied with Enzyme:

10X Antarctic Phosphatase Reaction Buffer.

**Reaction Conditions:** 1X Antarctic Phosphatase Reaction Buffer. Incubate at 37°C.

### 1X Antarctic Phosphatase Reaction Buffer:

50 mM Bis-Tris-Propane HCl

1 mM MgCl<sub>2</sub>

0.1 mM ZnCl<sub>2</sub>

pH 6.0 @ 25°C

**Unit Definition:** One unit is defined as the amount of enzyme that will dephosphorylate 1 µg of pUC19 vector DNA cut with HindIII (5' protruding ends), HincII (blunts ends) or Pst I (5' recessed ends) in 30 minutes at 37°C. Dephosphorylation is defined as > 95% inhibition of recirculation in a self-ligation reaction and is measured by transformation into *E. coli*.

**Unit Assay Conditions:** Vector DNA is dephosphorylated in restriction endonuclease buffer supplemented with Antarctic Phosphatase Reaction Buffer. Ligation is performed with 50 ng of vector using the NEB Quick Ligation Kit (NEB #M2200).

### Vector Dephosphorylation Protocol:

1. Add 1/10 volume of 10X Antarctic Phosphatase Reaction Buffer to 1–5 µg of DNA cut with any restriction endonuclease in any buffer.
2. Add 1 µl of Antarctic Phosphatase (5 units) and mix.
3. Incubate for 15 minutes at 37°C for 5' extensions or blunt-ends, 60 minutes for 3' extensions.
4. Heat inactivate for 5 minutes at 70°C (or as required to inactivate the restriction enzyme).
5. Proceed with ligation.

**Unit Definition:** One unit is defined as the amount of enzyme that will dephosphorylate 1 µg of pUC19 vector DNA cut with HindIII (5' protruding ends), HincII (blunts ends) or Pst I (5' recessed ends) in 30 minutes at 37°C. Dephosphorylation is defined as > 95% inhibition of recirculation in a self-ligation reaction and is measured by transformation into *E. coli*.

**Unit Assay Conditions:** Vector DNA is dephosphorylated in restriction endonuclease buffer supplemented with Antarctic Phosphatase Reaction Buffer. Ligation is performed with 50 ng of vector using the NEB Quick Ligation Kit (NEB #M2200).

### Vector Dephosphorylation Protocol:

1. Add 1/10 volume of 10X Antarctic Phosphatase Reaction Buffer to 1–5 µg of DNA cut with any restriction endonuclease in any buffer.
2. Add 1 µl of Antarctic Phosphatase (5 units) and mix.
3. Incubate for 15 minutes at 37°C for 5' extensions or blunt-ends, 60 minutes for 3' extensions.
4. Heat inactivate for 5 minutes at 70°C (or as required to inactivate the restriction enzyme).
5. Proceed with ligation.

**Usage Notes:** Antarctic Phosphatase is also active in NEBuffers 1, 2, 3 or 4 as well as the NEBuffer for EcoRI **ONLY** when supplemented with 10X Antarctic Phosphatase Reaction Buffer to a final concentration of 1X. Adding 1/10 volume of 10X Antarctic Phosphatase Reaction Buffer will provide the amount of Zn<sup>2+</sup> that the enzyme requires for activity.

### Quality Control Assays

**Exonuclease Activity:** Incubation of a 50 µl reaction containing 50 units of Antarctic Phosphatase with 1 µg of a mixture of single and double-stranded [<sup>3</sup>H] *E. coli* DNA for 4 hours at 37°C released < 0.1% of the total radioactivity.

**Endonuclease Activity:** Incubation of a 50 µl reaction containing 5 units of Antarctic Phosphatase with 1 µg of φX174 RF I DNA for 4 hours at 37°C resulted in < 10% conversion to RFI as determined by agarose gel electrophoresis.

**RNase Activity:** Incubation of a 10 µl reaction containing 50 units of Antarctic Phosphatase with 40 ng of fluorescein labeled RNA transcript for 4 hours at 37°C resulted in no detectable degradation of the RNA as determined by gel electrophoresis using fluorescence detection.

CERTIFICATE OF ANALYSIS

**Usage Notes:** Antarctic Phosphatase is also active in NEBuffers 1, 2, 3 or 4 as well as the NEBuffer for EcoRI **ONLY** when supplemented with 10X Antarctic Phosphatase Reaction Buffer to a final concentration of 1X. Adding 1/10 volume of 10X Antarctic Phosphatase Reaction Buffer will provide the amount of Zn<sup>2+</sup> that the enzyme requires for activity.

### Quality Control Assays

**Exonuclease Activity:** Incubation of a 50 µl reaction containing 50 units of Antarctic Phosphatase with 1 µg of a mixture of single and double-stranded [<sup>3</sup>H] *E. coli* DNA for 4 hours at 37°C released < 0.1% of the total radioactivity.

**Endonuclease Activity:** Incubation of a 50 µl reaction containing 5 units of Antarctic Phosphatase with 1 µg of φX174 RF I DNA for 4 hours at 37°C resulted in < 10% conversion to RFI as determined by agarose gel electrophoresis.

**RNase Activity:** Incubation of a 10 µl reaction containing 50 units of Antarctic Phosphatase with 40 ng of fluorescein labeled RNA transcript for 4 hours at 37°C resulted in no detectable degradation of the RNA as determined by gel electrophoresis using fluorescence detection.

CERTIFICATE OF ANALYSIS

**Physical Purity:** Purified to > 95% homogeneity as determined by SDS-PAGE analysis using Coomassie Blue detection.

**Transformation Assay:** pUC19 was cleaved with HindIII, HincII or PstI, each purified by Qiaprep™ spin column and resuspended in water at 0.2 mg/ml. 1 µg of each DNA was treated with 5 units of Antarctic Phosphatase. Each reaction was carried out with 1 µg of DNA in a 50 µl reaction volume with 1X Antarctic Phosphatase Reaction Buffer for 30 minutes at 37°C, followed by heat inactivation at 65°C for 5 minutes.

Ligations were performed using the NEB Quick Ligation Kit protocol with 2.5 µl (50 ng) of vector DNA directly from the heat-inactivated phosphatase reaction mix. Inserts were included in 3-fold molar excess. Either HindIII cleaved fragments of λ DNA, HaeIII cleaved fragments of φX174 RF I DNA or PstI cleaved fragments of λ DNA were inserted as indicated.

5 ng of each ligation was transformed into *E. coli* DH5-α. The equivalent of 1.0 ng were plated on LB plates that contained IPTG, X-gal and ampicillin.

Page 2 (M0289)

**Physical Purity:** Purified to > 95% homogeneity as determined by SDS-PAGE analysis using Coomassie Blue detection.

**Transformation Assay:** pUC19 was cleaved with HindIII, HincII or PstI, each purified by Qiaprep™ spin column and resuspended in water at 0.2 mg/ml. 1 µg of each DNA was treated with 5 units of Antarctic Phosphatase. Each reaction was carried out with 1 µg of DNA in a 50 µl reaction volume with 1X Antarctic Phosphatase Reaction Buffer for 30 minutes at 37°C, followed by heat inactivation at 65°C for 5 minutes.

Ligations were performed using the NEB Quick Ligation Kit protocol with 2.5 µl (50 ng) of vector DNA directly from the heat-inactivated phosphatase reaction mix. Inserts were included in 3-fold molar excess. Either HindIII cleaved fragments of λ DNA, HaeIII cleaved fragments of φX174 RF I DNA or PstI cleaved fragments of λ DNA were inserted as indicated.

5 ng of each ligation was transformed into *E. coli* DH5-α. The equivalent of 1.0 ng were plated on LB plates that contained IPTG, X-gal and ampicillin.

Page 2 (M0289)

|                                     | Antarctic Phosphatase Units | Blue Colonies | White Colonies | Background Reduction (%) | Insert (%) |
|-------------------------------------|-----------------------------|---------------|----------------|--------------------------|------------|
| <b>5' Overhang</b>                  |                             |               |                |                          |            |
| HindIII cleaved and ligated         | 0                           | 1360          | 0              |                          |            |
| Phosphatase treated and ligated     | 5                           | 16            | 0              |                          |            |
| Phosphatase treated, insert ligated | 5                           | 24            | 216            | 98.2                     | 15.9       |
| <b>Blunt End</b>                    |                             |               |                |                          |            |
| HincII cleaved and ligated          | 0                           | 1028          | 0              |                          |            |
| Phosphatase treated and ligated     | 5                           | 0             | 0              |                          |            |
| Phosphatase treated, insert ligated | 5                           | 24            | 92             | 97.7                     | 8.9        |
| <b>3' Overhang</b>                  |                             |               |                |                          |            |
| PstI cleaved and ligated            | 0                           | 872           | 0              |                          |            |
| Phosphatase treated and ligated     | 5                           | 8             | 0              |                          |            |
| Phosphatase treated, insert ligated | 5                           | 20            | 348            | 97.7                     | 40.0       |

**Heat Inactivation:** 70°C for 5 minutes.

**References:**

1. Sambrook, J., Fritsch, E.F. and Maniatis, T. (1989) *Molecular Cloning: A Laboratory Manual*, (2nd ed.), (p. 5.72). Cold Spring Harbor: Cold Spring Harbor Laboratory Press.
2. Rina, M. et al. (2000) *Eur. J. Biochem.* 267, 1230–1238.
3. Guthrie, E., unpublished observations

|                                     | Antarctic Phosphatase Units | Blue Colonies | White Colonies | Background Reduction (%) | Insert (%) |
|-------------------------------------|-----------------------------|---------------|----------------|--------------------------|------------|
| <b>5' Overhang</b>                  |                             |               |                |                          |            |
| HindIII cleaved and ligated         | 0                           | 1360          | 0              |                          |            |
| Phosphatase treated and ligated     | 5                           | 16            | 0              |                          |            |
| Phosphatase treated, insert ligated | 5                           | 24            | 216            | 98.2                     | 15.9       |
| <b>Blunt End</b>                    |                             |               |                |                          |            |
| HincII cleaved and ligated          | 0                           | 1028          | 0              |                          |            |
| Phosphatase treated and ligated     | 5                           | 0             | 0              |                          |            |
| Phosphatase treated, insert ligated | 5                           | 24            | 92             | 97.7                     | 8.9        |
| <b>3' Overhang</b>                  |                             |               |                |                          |            |
| PstI cleaved and ligated            | 0                           | 872           | 0              |                          |            |
| Phosphatase treated and ligated     | 5                           | 8             | 0              |                          |            |
| Phosphatase treated, insert ligated | 5                           | 20            | 348            | 97.7                     | 40.0       |

**Heat Inactivation:** 70°C for 5 minutes.

**References:**

1. Sambrook, J., Fritsch, E.F. and Maniatis, T. (1989) *Molecular Cloning: A Laboratory Manual*, (2nd ed.), (p. 5.72). Cold Spring Harbor: Cold Spring Harbor Laboratory Press.
2. Rina, M. et al. (2000) *Eur. J. Biochem.* 267, 1230–1238.
3. Guthrie, E., unpublished observations