

# RNase H



1-800-632-7799  
info@neb.com  
www.neb.com



M0297S 006170219021

## M0297S



250 units 5,000 U/ml Lot: 0061702

RECOMBINANT Store at -20°C Exp: 2/19

**Description:** Ribonuclease H (RNase) H is an endoribonuclease which specifically hydrolyzes the phosphodiester bonds of RNA which is hybridized to DNA. This enzyme does not digest single or double-stranded DNA.

**Source:** An *E. coli* strain that carries the cloned RNase H gene (rnh) from *Escherichia coli*

# RNase H



1-800-632-7799  
info@neb.com  
www.neb.com



M0297S 006170219021

## M0297S



250 units 5,000 U/ml Lot: 0061702

RECOMBINANT Store at -20°C Exp: 2/19

**Description:** Ribonuclease H (RNase) H is an endoribonuclease which specifically hydrolyzes the phosphodiester bonds of RNA which is hybridized to DNA. This enzyme does not digest single or double-stranded DNA.

**Source:** An *E. coli* strain that carries the cloned RNase H gene (rnh) from *Escherichia coli*

### Applications:

- Removal of poly(A) tails of mRNA hybridized to poly(dT)
- Removal of mRNA during second strand cDNA synthesis

Supplied in: 50 mM KCl, 10 mM Tris-HCl (pH 7.4 @ 25°C), 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol.

### Reagents Supplied with Enzyme:

10X RNase H Reaction Buffer.

**Reaction Conditions:** 1X RNase H Reaction Buffer. Incubate at 37°C.

### 1X RNase H Reaction Buffer:

75 mM KCl  
50 mM Tris-HCl  
3 mM MgCl<sub>2</sub>  
10 mM dithiothreitol  
pH 8.3 @ 25°C

### Applications:

- Removal of poly(A) tails of mRNA hybridized to poly(dT)
- Removal of mRNA during second strand cDNA synthesis

Supplied in: 50 mM KCl, 10 mM Tris-HCl (pH 7.4 @ 25°C), 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol.

### Reagents Supplied with Enzyme:

10X RNase H Reaction Buffer.

**Reaction Conditions:** 1X RNase H Reaction Buffer. Incubate at 37°C.

### 1X RNase H Reaction Buffer:

75 mM KCl  
50 mM Tris-HCl  
3 mM MgCl<sub>2</sub>  
10 mM dithiothreitol  
pH 8.3 @ 25°C

**Unit Definition:** One unit is defined as the amount of enzyme required to produce 1 nmol of ribonucleotides from 20 picomoles of a fluorescently labelled 50 base pair RNA-DNA hybrid in a total reaction volume of 50 µl in 20 minutes at 37°C.

### Quality Control Assays

**RNase Assay:** Incubation of a 10 µl reaction containing 50 units of RNase H with 40 ng of RNA transcripts for 1 hour at 37°C results in no detectable degradation of the RNA as determined by gel electrophoresis.

**SS DNA Exonuclease Activity:** Incubation of 50 units of enzyme with 1 µg of sonicated and denatured [<sup>3</sup>H]-DNA (10<sup>5</sup> cpm/ug) for 30 minutes at 37°C in 50 µl reaction buffer released < 0.1% radioactivity.

**Unit Definition:** One unit is defined as the amount of enzyme required to produce 1 nmol of ribonucleotides from 20 picomoles of a fluorescently labelled 50 base pair RNA-DNA hybrid in a total reaction volume of 50 µl in 20 minutes at 37°C.

### Quality Control Assays

**RNase Assay:** Incubation of a 10 µl reaction containing 50 units of RNase H with 40 ng of RNA transcripts for 1 hour at 37°C results in no detectable degradation of the RNA as determined by gel electrophoresis.

**SS DNA Exonuclease Activity:** Incubation of 50 units of enzyme with 1 µg of sonicated and denatured [<sup>3</sup>H]-DNA (10<sup>5</sup> cpm/ug) for 30 minutes at 37°C in 50 µl reaction buffer released < 0.1% radioactivity.

**Endonuclease Activity:** Incubation of a 50 µl reaction containing 50 units of RNase H with 1 µg of ϕX174 RF I DNA for 4 hours at 37°C resulted in < 10% conversion to RFI as determined by agarose gel electrophoresis.

**Heat Inactivation:** 65°C for 20 minutes.

### References:

1. Gubler, U. and Hoffman, B.J. (1983) *Gene* 25, 263–269.
2. Davis, R. et al. (1988) *Cell Biol.* 8, 4745–4755.
3. Donnis-Keller, H. (1979) *Nucleic Acid Res.* 7, 179.
4. Goodwin, E. C. and Rottman, F.M. (1992) *Nucleic Acids Res.* 20, 916.



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS

**Endonuclease Activity:** Incubation of a 50 µl reaction containing 50 units of RNase H with 1 µg of ϕX174 RF I DNA for 4 hours at 37°C resulted in < 10% conversion to RFI as determined by agarose gel electrophoresis.

**Heat Inactivation:** 65°C for 20 minutes.

### References:

1. Gubler, U. and Hoffman, B.J. (1983) *Gene* 25, 263–269.
2. Davis, R. et al. (1988) *Cell Biol.* 8, 4745–4755.
3. Donnis-Keller, H. (1979) *Nucleic Acid Res.* 7, 179.
4. Goodwin, E. C. and Rottman, F.M. (1992) *Nucleic Acids Res.* 20, 916.



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS