

## Shrimp Alkaline Phosphatase (rSAP)



1-800-632-7799  
info@neb.com  
www.neb.com



M0371S 001130615061

# M0371S



**500 units**      **1,000 U/ml**      **Lot: 0011306**  
**Store at -20°C**      **Exp: 6/15**

**Description:** Shrimp Alkaline Phosphatase (rSAP) is a heat labile alkaline phosphatase purified from a recombinant source. rSAP is identical to the native enzyme and contains no affinity tags or other modifications. rSAP nonspecifically catalyzes the dephosphorylation of 5' and 3' ends of DNA and RNA phosphomonoesters but does not degrade diphosphate and triphosphate linkages. Also, rSAP hydrolyses ribo-, as well as deoxyribonucleoside triphosphates (rNTPs and dNTPs). rSAP is useful in many molecular biology applications such as

## Shrimp Alkaline Phosphatase (rSAP)



1-800-632-7799  
info@neb.com  
www.neb.com



M0371S 001130615061

# M0371S



**500 units**      **1,000 U/ml**      **Lot: 0011306**  
**Store at -20°C**      **Exp: 6/15**

**Description:** Shrimp Alkaline Phosphatase (rSAP) is a heat labile alkaline phosphatase purified from a recombinant source. rSAP is identical to the native enzyme and contains no affinity tags or other modifications. rSAP nonspecifically catalyzes the dephosphorylation of 5' and 3' ends of DNA and RNA phosphomonoesters but does not degrade diphosphate and triphosphate linkages. Also, rSAP hydrolyses ribo-, as well as deoxyribonucleoside triphosphates (rNTPs and dNTPs). rSAP is useful in many molecular biology applications such as

the removal of phosphorylated ends of DNA and RNA for subsequent use in cloning or end-labeling of probes. In cloning, dephosphorylation prevents religation of linearized plasmid DNA. The enzyme acts on 5' protruding, 5' recessed and blunt ends. rSAP may also be used to degrade unincorporated dNTPs in PCR reactions to prepare templates for DNA sequencing or SNP analysis. rSAP is completely and irreversibly inactivated by heating at 65°C for 5 minutes, thereby making removal of rSAP prior to ligation or end-labeling unnecessary.

**Source:** A *Pichia pastoris* clone that carries the shrimp alkaline phosphatase gene from Northern shrimp *Pandalus borealis* (1,2).

**Specific Activity:** ~ 3,000 units/mg

**Molecular Weight:** rSAP is a homodimer. The molecular weight of the monomer is 54 kDa.

### Applications:

- Dephosphorylation of cloning vector DNA to prevent recircularization during ligation.
- Dephosphorylation of DNA prior to end-labeling using T4 Polynucleotide Kinase.
- Treatment of dNTPs in PCR reactions prior to

the removal of phosphorylated ends of DNA and RNA for subsequent use in cloning or end-labeling of probes. In cloning, dephosphorylation prevents religation of linearized plasmid DNA. The enzyme acts on 5' protruding, 5' recessed and blunt ends. rSAP may also be used to degrade unincorporated dNTPs in PCR reactions to prepare templates for DNA sequencing or SNP analysis. rSAP is completely and irreversibly inactivated by heating at 65°C for 5 minutes, thereby making removal of rSAP prior to ligation or end-labeling unnecessary.

**Source:** A *Pichia pastoris* clone that carries the shrimp alkaline phosphatase gene from Northern shrimp *Pandalus borealis* (1,2).

**Specific Activity:** ~ 3,000 units/mg

**Molecular Weight:** rSAP is a homodimer. The molecular weight of the monomer is 54 kDa.

### Applications:

- Dephosphorylation of cloning vector DNA to prevent recircularization during ligation.
- Dephosphorylation of DNA prior to end-labeling using T4 Polynucleotide Kinase.
- Treatment of dNTPs in PCR reactions prior to

sequencing or SNP analysis.

- Dephosphorylation of DNA and RNA.

Supplied in: 25 mM Tris-HCl (pH 7.5), 1 mM MgCl<sub>2</sub> and 50% glycerol.

### Reagents Supplied with Enzyme:

10X rSAP Reaction Buffer

### Reaction Conditions:

1X rSAP Reaction Buffer.

Incubate at 37°C.

### 1X rSAP Reaction Buffer:

20 mM Tris-HCl

10 mM MgCl<sub>2</sub>

pH 8.0 @ 25°C

**Unit Definition:** One unit is defined as the amount of enzyme that hydrolyzes 1 μmol of *p*-Nitrophenyl Phosphate, PNPP (NEB #P0757) in a total reaction volume of 1 ml in 1 minute at 37°C.

**Unit Assay Conditions:** 1 M diethanolamine-HCl (pH 9.8), 0.5 mM MgCl<sub>2</sub>, 50 mM *p*-Nitrophenyl Phosphate. These conditions are only used for quantitating enzyme activity.

sequencing or SNP analysis.

- Dephosphorylation of DNA and RNA.

Supplied in: 25 mM Tris-HCl (pH 7.5), 1 mM MgCl<sub>2</sub> and 50% glycerol.

### Reagents Supplied with Enzyme:

10X rSAP Reaction Buffer

### Reaction Conditions:

1X rSAP Reaction Buffer.

Incubate at 37°C.

### 1X rSAP Reaction Buffer:

20 mM Tris-HCl

10 mM MgCl<sub>2</sub>

pH 8.0 @ 25°C

**Unit Definition:** One unit is defined as the amount of enzyme that hydrolyzes 1 μmol of *p*-Nitrophenyl Phosphate, PNPP (NEB #P0757) in a total reaction volume of 1 ml in 1 minute at 37°C.

**Unit Assay Conditions:** 1 M diethanolamine-HCl (pH 9.8), 0.5 mM MgCl<sub>2</sub>, 50 mM *p*-Nitrophenyl Phosphate. These conditions are only used for quantitating enzyme activity.

**Functional Assay:** Dephosphorylation with rSAP of a restriction enzyme-digested vector DNA with 5' recessed ends, the least favorable type for dephosphorylation, reduces re-ligation to < 0.5% compared to untreated control as measured by transformation into *E. coli*. rSAP has been functionally tested in the following protocol:

### Protocol for Dephosphorylation of 5'-ends of DNA

1. Resuspend 1 pmol of DNA ends (about 1 μg of a 3 kb plasmid) in purified water.
2. Prepare a 20 μl reaction as follows:

DNA	≥ 1 μl
rSAP Reaction Buffer (10X)	2 μl
rSAP (1 unit/μl)	1 μl
H <sub>2</sub> O, purified	to 20 μl

*Note: Scale larger reaction volumes proportionally.*
3. Incubate at 37°C for 30 minutes.
4. Stop reaction by heat-inactivation at 65°C for 5 minutes.

(see other side)

CERTIFICATE OF ANALYSIS

**Functional Assay:** Dephosphorylation with rSAP of a restriction enzyme-digested vector DNA with 5' recessed ends, the least favorable type for dephosphorylation, reduces re-ligation to < 0.5% compared to untreated control as measured by transformation into *E. coli*. rSAP has been functionally tested in the following protocol:

### Protocol for Dephosphorylation of 5'-ends of DNA

1. Resuspend 1 pmol of DNA ends (about 1 μg of a 3 kb plasmid) in purified water.
2. Prepare a 20 μl reaction as follows:

DNA	≥ 1 μl
rSAP Reaction Buffer (10X)	2 μl
rSAP (1 unit/μl)	1 μl
H <sub>2</sub> O, purified	to 20 μl

*Note: Scale larger reaction volumes proportionally.*
3. Incubate at 37°C for 30 minutes.
4. Stop reaction by heat-inactivation at 65°C for 5 minutes.

(see other side)

CERTIFICATE OF ANALYSIS

### Protocol for Dephosphorylation of 5'-ends of DNA in Restriction Enzyme Reaction

- Digest 1–5 µg of plasmid DNA in a 20 µl reaction as follows:

DNA	≥ 1 µl
Restriction Enzyme Buffer (10X)	2 µl
Restriction Endonuclease	1 µl
H <sub>2</sub> O, purified	to 20 µl

*Note: Scale larger reaction volumes proportionally.*
- Incubate at 37°C for 60 minutes or follow manufacturer's recommendations.
- Add 1 unit of rSAP for every 1 pmol of DNA ends (about 1 µg of a 3 kb plasmid) and incubate at 37°C for 30–60 minutes.
- Stop reaction by heat-inactivation of rSAP and restriction enzyme (follow manufacturer's recommendations).
- Proceed with ligation.

*Note: If restriction enzyme cannot be heat-inactivated, DNA purification is required before ligation.*

Page 2 (M0371)

### Protocol for Dephosphorylation of 5'-ends of DNA in Restriction Enzyme Reaction

- Digest 1–5 µg of plasmid DNA in a 20 µl reaction as follows:

DNA	≥ 1 µl
Restriction Enzyme Buffer (10X)	2 µl
Restriction Endonuclease	1 µl
H <sub>2</sub> O, purified	to 20 µl

*Note: Scale larger reaction volumes proportionally.*
- Incubate at 37°C for 60 minutes or follow manufacturer's recommendations.
- Add 1 unit of rSAP for every 1 pmol of DNA ends (about 1 µg of a 3 kb plasmid) and incubate at 37°C for 30–60 minutes.
- Stop reaction by heat-inactivation of rSAP and restriction enzyme (follow manufacturer's recommendations).
- Proceed with ligation.

*Note: If restriction enzyme cannot be heat-inactivated, DNA purification is required before ligation.*

Page 2 (M0371)

### Usage Notes:

- rSAP, as are most alkaline phosphatases, is a Zn<sup>2+</sup> and Mg<sup>2+</sup>-dependent enzyme. Our formulation has tightly bound zinc atoms in the active center and does not require supplemental zinc or other additives.
- rSAP is also active in 1X NEBuffers 1.1, 2.1, 3.1 and CutSmart™ NEBuffer.
- rSAP activity is enhanced in the presence of monovalent salts.
- rSAP is inhibited by metal chelators (e.g. EDTA), inorganic phosphate and phosphate analogs.
- The rSAP activity is decreased in the presence of reducing agents (DTT, β-ME).

### Quality Controls Assays

**Exonuclease Activity:** Incubation of a 50 µl reaction containing 10 units of rSAP with 1 µg of a mixture of single and double-stranded [<sup>3</sup>H] *E. coli* DNA for 4 hours at 37°C released < 0.1% of the total radioactivity.

### Usage Notes:

- rSAP, as are most alkaline phosphatases, is a Zn<sup>2+</sup> and Mg<sup>2+</sup>-dependent enzyme. Our formulation has tightly bound zinc atoms in the active center and does not require supplemental zinc or other additives.
- rSAP is also active in 1X NEBuffers 1.1, 2.1, 3.1 and CutSmart™ NEBuffer.
- rSAP activity is enhanced in the presence of monovalent salts.
- rSAP is inhibited by metal chelators (e.g. EDTA), inorganic phosphate and phosphate analogs.
- The rSAP activity is decreased in the presence of reducing agents (DTT, β-ME).

### Quality Controls Assays

**Exonuclease Activity:** Incubation of a 50 µl reaction containing 10 units of rSAP with 1 µg of a mixture of single and double-stranded [<sup>3</sup>H] *E. coli* DNA for 4 hours at 37°C released < 0.1% of the total radioactivity.

**Endonuclease Activity:** Incubation of a 50 µl reaction containing 5 units of rSAP with 1 µg of φX174 RF I DNA for 4 hours at 37°C resulted in < 10% conversion to RF II as determined by agarose gel electrophoresis.

**RNase Activity:** Incubation of a 10 µl reaction containing 5 units of rSAP with 40 ng of fluorescein labeled RNA transcript for 4 hours at 37°C resulted in < 10% degradation of the RNA as determined by gel electrophoresis using fluorescence detection.

**Physical Purity:** Purified to > 95% homogeneity as determined by SDS-PAGE analysis using Coomassie Blue staining detection.

**Heat Inactivation:** 65°C for 5 minutes

### References:

- Olsen, R. L. et al. (1991) *Comp. Biochem. Physiol.* 99B(4): 755–761.
- Nilsen, I.W. et al. (2001) *Comp. Biochem. Physiol.* 129B(4): 853–861.

**Endonuclease Activity:** Incubation of a 50 µl reaction containing 5 units of rSAP with 1 µg of φX174 RF I DNA for 4 hours at 37°C resulted in < 10% conversion to RF II as determined by agarose gel electrophoresis.

**RNase Activity:** Incubation of a 10 µl reaction containing 5 units of rSAP with 40 ng of fluorescein labeled RNA transcript for 4 hours at 37°C resulted in < 10% degradation of the RNA as determined by gel electrophoresis using fluorescence detection.

**Physical Purity:** Purified to > 95% homogeneity as determined by SDS-PAGE analysis using Coomassie Blue staining detection.

**Heat Inactivation:** 65°C for 5 minutes

### References:

- Olsen, R. L. et al. (1991) *Comp. Biochem. Physiol.* 99B(4): 755–761.
- Nilsen, I.W. et al. (2001) *Comp. Biochem. Physiol.* 129B(4): 853–861.

### Companion Products Sold Separately:

T4 DNA Ligase	
#M0202S	20,000 units
#M0202L	100,000 units
#M0202T	20,000 units
#M0202M	100,000 units

Quick Ligation™ Kit	
#M2200S	30 rxns
#M2200L	150 rxns

Instant Sticky-end Ligase Master Mix	
#M0370S	50 rxns
#M0370L	250 rxns

Blunt/TA Ligase Master Mix	
#M0367S	50 rxns
#M0367L	250 rxns



New England Biolabs® is a registered trademark of New England Biolabs, Inc.

CUTSMART™ and QUICK LIGATION™ are trademarks of New England Biolabs, Inc.

rSAP is sourced from ArcticZymes AS; Quality assured, packaged and distributed by New England Biolabs, Inc.

### Companion Products Sold Separately:

T4 DNA Ligase	
#M0202S	20,000 units
#M0202L	100,000 units
#M0202T	20,000 units
#M0202M	100,000 units

Quick Ligation™ Kit	
#M2200S	30 rxns
#M2200L	150 rxns

Instant Sticky-end Ligase Master Mix	
#M0370S	50 rxns
#M0370L	250 rxns

Blunt/TA Ligase Master Mix	
#M0367S	50 rxns
#M0367L	250 rxns



New England Biolabs® is a registered trademark of New England Biolabs, Inc.

CUTSMART™ and QUICK LIGATION™ are trademarks of New England Biolabs, Inc.

rSAP is sourced from ArcticZymes AS; Quality assured, packaged and distributed by New England Biolabs, Inc.