

Oligo d(T)<sub>23</sub> VN



1-800-632-7799  
info@neb.com  
www.neb.com



S1327S 010151018101

**S1327S**

**1.0 A<sub>260</sub> unit**      **Lot: 0101510**  
**Store at -20°C**      **Exp: 10/18**

**Description:** Oligo d(T)<sub>23</sub> VN is used for the priming and sequencing of mRNA adjacent to the 3' -poly A tail (1).

V = A or G or C  
N = A or G or C or T

Supplied as a lyophilized triethylammonium salt.  
Reconstitute with water.

**This Primer Does Not Contain a 5' Phosphate.**

Oligo d(T)<sub>23</sub> VN



1-800-632-7799  
info@neb.com  
www.neb.com



S1327S 010151018101

**S1327S**

**1.0 A<sub>260</sub> unit**      **Lot: 0101510**  
**Store at -20°C**      **Exp: 10/18**

**Description:** Oligo d(T)<sub>23</sub> VN is used for the priming and sequencing of mRNA adjacent to the 3' -poly A tail (1).

V = A or G or C  
N = A or G or C or T

Supplied as a lyophilized triethylammonium salt.  
Reconstitute with water.

**This Primer Does Not Contain a 5' Phosphate.**

**Chromatographic Analysis:** HPLC C8 Reverse Phse: 10 µM, 10cm x 8 mm Waters Z-Module Rad-Pak, 0.1M TEAB, (pH7.0) to 15% acetonitrile, 45 minute linear gradient. Elution time 30.5 minutes.

**Enzymatic Analysis:** Results of electrophoretogram (20% acrylamide gel)

Primer <sup>32</sup>P labelled with T4 Polynucleotide Kinase showing one homogeneous band.

**Special Data** = Lambda max 258 nm.

µg/OD = 37.6

**1.0 A<sub>260</sub> unit** = 4.97 nmoles  
Dissolve 99.4 µl dH<sub>2</sub>O to make a 50 µM solution.

**References:**

1. Khan S. A. et al. (1991) *Nucleic Acids Res.* 19, 1715
2. Liao, J. and Gong, Z. (1997) *Biotechniques* 23, 368-370.

**Chromatographic Analysis:** HPLC C8 Reverse Phse: 10 µM, 10cm x 8 mm Waters Z-Module Rad-Pak, 0.1M TEAB, (pH7.0) to 15% acetonitrile, 45 minute linear gradient. Elution time 30.5 minutes.

**Enzymatic Analysis:** Results of electrophoretogram (20% acrylamide gel)

Primer <sup>32</sup>P labelled with T4 Polynucleotide Kinase showing one homogeneous band.

**Special Data** = Lambda max 258 nm.

µg/OD = 37.6

**1.0 A<sub>260</sub> unit** = 4.97 nmoles  
Dissolve 99.4 µl dH<sub>2</sub>O to make a 50 µM solution.

**References:**

1. Khan S. A. et al. (1991) *Nucleic Acids Res.* 19, 1715
2. Liao, J. and Gong, Z. (1997) *Biotechniques* 23, 368-370.



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS