
New England Biolabs Product Specification

Product Name: Monarch[®] DNA Wash Buffer
Catalog #: T1032L
Shelf Life: 24 months
Storage Temp: 25°C
Composition: Proprietary
Specification Version: PS-T1032L v1.0
Effective Date: 03 Dec 2017

Assay Name/Specification (minimum release criteria)

Functional Testing (Fragment Purification and Restriction Digestion) - A linear pLA3 fragment is added to 8 Monarch [®] DNA Cleanup Columns, eluted and digested with DraIII-HF [®] resulting in ≥90% conversion to two-fragment pattern observed by gel electrophoresis in ≥75% of the samples.



Date 03 Dec 2017

Derek Robinson
Director of Quality Control

