

---

## New England Biolabs Product Specification

*Product Name:* Monarch<sup>®</sup> 2 ml Tubes  
*Catalog #:* T3003L  
*Storage Temp:* 25°C  
*Specification Version:* PS-T3003L v1.0  
*Effective Date:* 16 Dec 2020

<b>Assay Name/Specification (minimum release criteria)</b>
--

<b>Integrity (Physical Inspection)</b> - A representative selection of each production batch of Monarch <sup>®</sup> 2 ml Tubes is visually inspected for signs of damage or manufacturing defect prior to batch release.
---

One or more products referenced in this document may be covered by a 3rd-party trademark.  
Please visit [www.neb.com/trademarks](http://www.neb.com/trademarks) for additional information.



Date 16 Dec 2020

---

Derek Robinson  
Director, Quality Control

