

ThermoPol® Reaction Buffer Pack



1-800-632-7799
info@neb.com
www.neb.com



B9004S 001141119111

B9004S

6.0 ml (4 x 1.5 ml)

Lot: 0011411

Store at -20°C

Exp: 11/19

Description: 10X ThermoPol Reaction Buffer is optimized for use with Vent_R® and Deep Vent_R™ DNA Polymerases. This buffer also provides superior reaction conditions for other thermophilic DNA polymerases, including *Taq* DNA Polymerase, as well as various other DNA and RNA modifying enzymes. Supplemental MgSO₄ is provided to maximize flexibility for users that require additional Mg²⁺ for their specific application. See individual product descriptions for the recommended buffer for each enzyme, or visit the FAQ page for

this product on the NEB website for a complete list of enzymes for which this buffer is either recommended or compatible.

Reagents Supplied:

10X ThermoPol Reaction Buffer 4 x 1.5 ml
100 mM MgSO₄ 1.5 ml

Reaction Conditions: At a 1X concentration, the reaction buffer assures optimal activity of Vent_R, Deep Vent_R, *Taq* and other DNA polymerases, as well as other enzymes for which ThermoPol Reaction Buffer is recommended.

1X ThermoPol Reaction Buffer:

20 mM Tris-HCl
10 mM (NH₄)₂SO₄
10 mM KCl
2 mM MgSO₄
0.1% Triton® X-100
pH 8.8 @ 25°C

Quality Assurance

This reaction buffer and supplement are free of detectable nucleases.



NEW ENGLAND BIOLABS®, VENT® and THERMOPOL® are registered trademarks of New England Biolabs, Inc.

DEEP VENT™ is a trademark of New England Biolabs, Inc.

TRITON® is a registered trademark of Union Carbide Corporation.

CERTIFICATE OF ANALYSIS

ThermoPol® Reaction Buffer Pack



1-800-632-7799
info@neb.com
www.neb.com



B9004S 001141119111

B9004S

6.0 ml (4 x 1.5 ml)

Lot: 0011411

Store at -20°C

Exp: 11/19

Description: 10X ThermoPol Reaction Buffer is optimized for use with Vent_R® and Deep Vent_R™ DNA Polymerases. This buffer also provides superior reaction conditions for other thermophilic DNA polymerases, including *Taq* DNA Polymerase, as well as various other DNA and RNA modifying enzymes. Supplemental MgSO₄ is provided to maximize flexibility for users that require additional Mg²⁺ for their specific application. See individual product descriptions for the recommended buffer for each enzyme, or visit the FAQ page for

this product on the NEB website for a complete list of enzymes for which this buffer is either recommended or compatible.

Reagents Supplied:

10X ThermoPol Reaction Buffer 4 x 1.5 ml
100 mM MgSO₄ 1.5 ml

Reaction Conditions: At a 1X concentration, the reaction buffer assures optimal activity of Vent_R, Deep Vent_R, *Taq* and other DNA polymerases, as well as other enzymes for which ThermoPol Reaction Buffer is recommended.

1X ThermoPol Reaction Buffer:

20 mM Tris-HCl
10 mM (NH₄)₂SO₄
10 mM KCl
2 mM MgSO₄
0.1% Triton® X-100
pH 8.8 @ 25°C

Quality Assurance

This reaction buffer and supplement are free of detectable nucleases.



NEW ENGLAND BIOLABS®, VENT® and THERMOPOL® are registered trademarks of New England Biolabs, Inc.

DEEP VENT™ is a trademark of New England Biolabs, Inc.

TRITON® is a registered trademark of Union Carbide Corporation.

CERTIFICATE OF ANALYSIS