

AluI Methyltransferase



1-800-632-7799
info@neb.com
www.neb.com




M0220S 008160718071

M0220S



100 units 5,000 U/ml Lot: 0081607
RECOMBINANT Store at -20°C Exp: 7/18

Methylation Site:



Description: AluI Methyltransferase modifies the cytosine residue (C⁵) in the sequence above.

Source: An *E. coli* strain that carries the cloned AluI modification gene from *Arthrobacter luteus* (ATCC 21606).

AluI Methyltransferase



1-800-632-7799
info@neb.com
www.neb.com




M0220S 008160718071

M0220S



100 units 5,000 U/ml Lot: 0081607
RECOMBINANT Store at -20°C Exp: 7/18

Methylation Site:



Description: AluI Methyltransferase modifies the cytosine residue (C⁵) in the sequence above.

Source: An *E. coli* strain that carries the cloned AluI modification gene from *Arthrobacter luteus* (ATCC 21606).

Supplied in: 50 mM KCl, 10 mM Tris-HCl (pH 7.5), 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol.

Reagents Supplied with Enzyme:

10X AluI Methyltransferase Reaction Buffer,
400X S-adenosylmethionine (32 mM).

Reaction Conditions: 1X AluI Methyltransferase Reaction Buffer, 80 µM S-adenosylmethionine.

1X AluI Methyltransferase Reaction Buffer:

50 mM Tris-HCl
10 mM EDTA
5 mM 2-mercaptoethanol
pH 7.5 @ 25°C

Protection Assay Conditions: AluI Methyltransferase is incubated with 1 µg of λ DNA in 10 µl 1X AluI Methyltransferase Reaction Buffer, supplemented with 80 µM S-adenosylmethionine, for one hour at 37°C followed by 15 minutes at 65°C. The extent of protection by AluI Methyltransferase is determined by the addition of 40 µl NEBuffer 1 supplemented with 10 mM MgCl₂ and 10 units of AluI restriction endonuclease. Incubation at 37°C for 30 minutes is followed by analysis on an agarose gel.

Supplied in: 50 mM KCl, 10 mM Tris-HCl (pH 7.5), 0.1 mM EDTA, 1 mM dithiothreitol, 200 µg/ml BSA and 50% glycerol.

Reagents Supplied with Enzyme:

10X AluI Methyltransferase Reaction Buffer,
400X S-adenosylmethionine (32 mM).

Reaction Conditions: 1X AluI Methyltransferase Reaction Buffer, 80 µM S-adenosylmethionine.

1X AluI Methyltransferase Reaction Buffer:

50 mM Tris-HCl
10 mM EDTA
5 mM 2-mercaptoethanol
pH 7.5 @ 25°C

Protection Assay Conditions: AluI Methyltransferase is incubated with 1 µg of λ DNA in 10 µl 1X AluI Methyltransferase Reaction Buffer, supplemented with 80 µM S-adenosylmethionine, for one hour at 37°C followed by 15 minutes at 65°C. The extent of protection by AluI Methyltransferase is determined by the addition of 40 µl NEBuffer 1 supplemented with 10 mM MgCl₂ and 10 units of AluI restriction endonuclease. Incubation at 37°C for 30 minutes is followed by analysis on an agarose gel.

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 µg of λ DNA in 1 hour at 37°C in a total reaction volume of 10 µl against cleavage by AluI restriction endonuclease.

Quality Control Assays

16-Hour Incubation: Incubation of 40 units of Methyltransferase with 1 µg λ DNA in 50 µl 1X NEBuffer 2 at 37°C resulted in no detectable endonuclease contamination.

Exonuclease Activity: Incubation of 40 units of AluI Methyltransferase with 1 µg sonicated ³H DNA (10⁵ cpm/µg) for 4 hours at 37°C in 50 µl NEBuffer 2 [50 mM NaCl, 10 mM Tris-HCl (pH 7.9 @ 25°C), 10 mM MgCl₂, 1 mM dithiothreitol] released 0.2% of the total radioactivity.

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 µg of λ DNA in 1 hour at 37°C in a total reaction volume of 10 µl against cleavage by AluI restriction endonuclease.

Quality Control Assays

16-Hour Incubation: Incubation of 40 units of Methyltransferase with 1 µg λ DNA in 50 µl 1X NEBuffer 2 at 37°C resulted in no detectable endonuclease contamination.

Exonuclease Activity: Incubation of 40 units of AluI Methyltransferase with 1 µg sonicated ³H DNA (10⁵ cpm/µg) for 4 hours at 37°C in 50 µl NEBuffer 2 [50 mM NaCl, 10 mM Tris-HCl (pH 7.9 @ 25°C), 10 mM MgCl₂, 1 mM dithiothreitol] released 0.2% of the total radioactivity.

Storage of SAM: S-adenosylmethionine is stored at -20°C as a 32 mM solution dissolved in 0.005 M sulfuric acid and 10% ethanol. Under these conditions SAM is stable for up to 6 months. SAM is unstable at (pH 7.5), 37°C, (1) and should be replenished in reactions incubated longer than 4 hours.

Methylation can be optimized by using fresh SAM.

Reference:

- Hoffman, J. L. (1986) *Biochemistry* 25, 4444-4449.

Companion Product:

S-adenosylmethionine (SAM)
#B9003S 0.5 ml



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS

Storage of SAM: S-adenosylmethionine is stored at -20°C as a 32 mM solution dissolved in 0.005 M sulfuric acid and 10% ethanol. Under these conditions SAM is stable for up to 6 months. SAM is unstable at (pH 7.5), 37°C, (1) and should be replenished in reactions incubated longer than 4 hours.

Methylation can be optimized by using fresh SAM.

Reference:

- Hoffman, J. L. (1986) *Biochemistry* 25, 4444-4449.

Companion Product:

S-adenosylmethionine (SAM)
#B9003S 0.5 ml



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS