

α 1-2 Fucosidase



1-800-632-7799
info@neb.com
www.neb.com



P0724S 011150816081

P0724S

1,000 units 20,000 U/ml Lot: 0111508
Store at 4°C Exp: 8/16

Description: α 1-2 Fucosidase is a highly specific exoglycosidase that catalyzes the hydrolysis of linear α 1-2 linked L-fucopyranosyl residues from oligosaccharides (1). In this case, a linear substrate is defined as having no branching on the adjacent residue.

New Reaction Buffer

α 1-2 Fucosidase



1-800-632-7799
info@neb.com
www.neb.com



P0724S 011150816081


P0724S

1,000 units 20,000 U/ml Lot: 0111508
Store at 4°C Exp: 8/16

Description: α 1-2 Fucosidase is a highly specific exoglycosidase that catalyzes the hydrolysis of linear α 1-2 linked L-fucopyranosyl residues from oligosaccharides (1). In this case, a linear substrate is defined as having no branching on the adjacent residue.

New Reaction Buffer


Specificity:




Note: p-nitrophenyl- α -L-fucopyranoside is NOT a substrate for this enzyme.

Detailed Specificity: Specificity can vary depending on incubation time and concentration of substrate (Fig 1).


A. 0.1 nm/ μ l substrate, 1 hour incubation




B. 0.1 nm/ μ l substrate, 1 hour incubation



C. 0.05 nm/ μ l substrate, 1 hour incubation




Specificity:




Note: p-nitrophenyl- α -L-fucopyranoside is NOT a substrate for this enzyme.

Detailed Specificity: Specificity can vary depending on incubation time and concentration of substrate (Fig 1).


A. 0.1 nm/ μ l substrate, 1 hour incubation



B. 0.1 nm/ μ l substrate, 1 hour incubation



C. 0.05 nm/ μ l substrate, 1 hour incubation



D. 0.1 nm/ μ l substrate, 1 hour incubation




Figure 1: Detailed specificity of α 1-2 Fucosidase. Reactions (A), (B) and (D) contained 1X G4 reaction buffer and reactions (C) contained 1X G6 reaction buffer. All reactions contain 1X BSA in a total reaction volume of 10 μ l, and all reactions were incubated at 37°C. All reactions contained 20 units of α 1-2 Fucosidase. Reaction (C) also contained 20 units of α -N-Acetylgalactosaminidase (#P0734). In reaction (C), the branched α 1-2 fucose is removed in the presence of both enzymes, but not by α 1-2 Fucosidase alone.

Source: Cloned from *Xanthomonas manihotis* and expressed in *E. coli*.

Supplied in: 50 mM NaCl, 20 mM Tris-HCl (pH 7.5 @ 25°C) and 0.1 mM Na₂EDTA.

D. 0.1 nm/ μ l substrate, 1 hour incubation




Figure 1: Detailed specificity of α 1-2 Fucosidase. Reactions (A), (B) and (D) contained 1X G4 reaction buffer and reactions (C) contained 1X G6 reaction buffer. All reactions contain 1X BSA in a total reaction volume of 10 μ l, and all reactions were incubated at 37°C. All reactions contained 20 units of α 1-2 Fucosidase. Reaction (C) also contained 20 units of α -N-Acetylgalactosaminidase (#P0734). In reaction (C), the branched α 1-2 fucose is removed in the presence of both enzymes, but not by α 1-2 Fucosidase alone.

Source: Cloned from *Xanthomonas manihotis* and expressed in *E. coli*.

Supplied in: 50 mM NaCl, 20 mM Tris-HCl (pH 7.5 @ 25°C) and 0.1 mM Na₂EDTA.

Reagents Supplied with Enzyme:

10X GlycoBuffer 1
100X BSA

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 μ g/ml BSA. Incubate at 37°C.

Optimal incubation times and enzyme concentrations must be determined empirically for a particular substrate.

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the α -L-fucose from 1 nmol of Fuc α 1-2Gal β 1-4Glc-7-amino-4-methyl-coumarin (AMC), in 1 hour at 37°C in a total reaction volume of 10 μ l.

Specific Activity: ~ 2,000,000 units/mg

Molecular Weight: 70,000 daltons

Unit Definition Assay: Two fold dilutions of α 1-2 Fucosidase are incubated with 1 nmol AMC-labeled substrate in 1X GlycoBuffer 1, supplemented with 100 μ g/ml BSA, in a 10 μ l reaction. The reaction mix is incubated for 1 hour at 37°C. Separation of reaction products are visualized via thin layer chromatography (1). (see other side)

CERTIFICATE OF ANALYSIS

Reagents Supplied with Enzyme:

10X GlycoBuffer 1
100X BSA

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 μ g/ml BSA. Incubate at 37°C.

Optimal incubation times and enzyme concentrations must be determined empirically for a particular substrate.

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the α -L-fucose from 1 nmol of Fuc α 1-2Gal β 1-4Glc-7-amino-4-methyl-coumarin (AMC), in 1 hour at 37°C in a total reaction volume of 10 μ l.

Specific Activity: ~ 2,000,000 units/mg

Molecular Weight: 70,000 daltons

Unit Definition Assay: Two fold dilutions of α 1-2 Fucosidase are incubated with 1 nmol AMC-labeled substrate in 1X GlycoBuffer 1, supplemented with 100 μ g/ml BSA, in a 10 μ l reaction. The reaction mix is incubated for 1 hour at 37°C. Separation of reaction products are visualized via thin layer chromatography (1). (see other side)

CERTIFICATE OF ANALYSIS

Quality Assurance: No contaminating exoglycosidase or proteolytic activity could be detected.

Quality Controls

Glycosidase Assays: 100 units of α 1-2 Fucosidase were incubated with 0.1 mM of fluorescently-labeled oligosaccharides and glycopeptides, in a 10 μ l reaction for 20 hours at 37°C. The reaction products were analyzed by TLC for digestion of substrate.

Physical Purity: Purified to > 95% homogeneity as determined by SDS-PAGE analysis using Coomassie Blue detection.

No other glycosidase activities were detected (ND) with the following substrates:

β -N-Acetyl-glucosaminidase:
GlcNAc β 1-4GlcNAc β 1-4GlcNAc-AMC ND

α -Fucosidase:
Gal β 1-4(Fuc α 1-3)GlcNAc β 1-3Gal β 1-4Glc-AMC ND

Page 2 (P0724)

Quality Assurance: No contaminating exoglycosidase or proteolytic activity could be detected.

Quality Controls

Glycosidase Assays: 100 units of α 1-2 Fucosidase were incubated with 0.1 mM of fluorescently-labeled oligosaccharides and glycopeptides, in a 10 μ l reaction for 20 hours at 37°C. The reaction products were analyzed by TLC for digestion of substrate.

Physical Purity: Purified to > 95% homogeneity as determined by SDS-PAGE analysis using Coomassie Blue detection.

No other glycosidase activities were detected (ND) with the following substrates:

β -N-Acetyl-glucosaminidase:
GlcNAc β 1-4GlcNAc β 1-4GlcNAc-AMC ND

α -Fucosidase:
Gal β 1-4(Fuc α 1-3)GlcNAc β 1-3Gal β 1-4Glc-AMC ND

Page 2 (P0724)

β -Galactosidase:
Gal β 1-3GlcNAc β 1-4Gal β 1-4Glc-AMC ND

α -Galactosidase:
Gal α 1-3Gal β 1-4GlcNAc-AMC ND

α -Neuraminidase:
Neu5Ac α 2-3Gal β 1-3GlcNAc β 1-3Gal β 1-4Glc-AMC ND

α -Mannosidase:
Man α 1-3Man β 1-4GlcNAc-AMC
Man α 1-6Man α 1-6(Man α 1-3)Man-AMC ND

β -Glucosidase:
Glc β 1-4Glc β 1-4Glc-AMC ND

β -Xylosidase:
Xyl β 1-4Xyl β 1-4Xyl β 1-4Xyl-AMC ND

β -Mannosidase:
Man β 1-4Man β 1-4Man-AMC ND

β -Galactosidase:
Gal β 1-3GlcNAc β 1-4Gal β 1-4Glc-AMC ND

α -Galactosidase:
Gal α 1-3Gal β 1-4GlcNAc-AMC ND

α -Neuraminidase:
Neu5Ac α 2-3Gal β 1-3GlcNAc β 1-3Gal β 1-4Glc-AMC ND

α -Mannosidase:
Man α 1-3Man β 1-4GlcNAc-AMC
Man α 1-6Man α 1-6(Man α 1-3)Man-AMC ND

β -Glucosidase:
Glc β 1-4Glc β 1-4Glc-AMC ND

β -Xylosidase:
Xyl β 1-4Xyl β 1-4Xyl β 1-4Xyl-AMC ND

β -Mannosidase:
Man β 1-4Man β 1-4Man-AMC ND

Endo F₁, F₂, H:
Dansylated invertase high mannose. ND

Endo F₂, F₃:
Dansylated fibrinogen biantennary. ND

PNGase F:
Fluoresceinated fetuin triantennary. ND

Protease Assay: After incubation of 100 units of α 1-2 Fucosidase with 0.2 nmol of a standard mixture of proteins in a 20 μ l reaction, for 20 hours at 37°C, no proteolytic activity could be detected by SDS-PAGE.

Note: Repeated freeze/thaw cycles may reduce activity. Recommended storage temperature has changed to 4°C.

Endo F₁, F₂, H:
Dansylated invertase high mannose. ND

Endo F₂, F₃:
Dansylated fibrinogen biantennary. ND

PNGase F:
Fluoresceinated fetuin triantennary. ND

Protease Assay: After incubation of 100 units of α 1-2 Fucosidase with 0.2 nmol of a standard mixture of proteins in a 20 μ l reaction, for 20 hours at 37°C, no proteolytic activity could be detected by SDS-PAGE.

Note: Repeated freeze/thaw cycles may reduce activity. Recommended storage temperature has changed to 4°C.

Reference:

1. Wong-Madden, S.T. and Landry, D. (1995) *Glycobiology* 5, 19–28.



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

Reference:

1. Wong-Madden, S.T. and Landry, D. (1995) *Glycobiology* 5, 19–28.



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.