

New England Biolabs Product Specification

Product Name: Monarch[®] PCR & DNA Cleanup Kit (5 µg)
Catalog #: T1030S/L
Kit Components: Monarch[®] DNA Elution Buffer (T1016)
Monarch[®] DNA Cleanup Binding Buffer (T1031)
Monarch[®] DNA Wash Buffer (T1032)
Monarch[®] DNA Cleanup Columns (5 µg) (T1034)
Monarch[®] Collection Tubes (2 ml) (T1018)
Shelf Life: 18 months
Storage Temp: 25°C
Specification Version: PS-T1030S/L v1.0
Effective Date: 03 Dec 2017

Assay Name/Specification (minimum release criteria)

Functional Testing (DNA Recovery and Purity) - Twenty-four Monarch[®] DNA Cleanup Columns are tested with 1 µg of 1 Kb Ladder resulting in ≥75% recovery in ≥80% of the samples. OD 260/280 and 260/230 are ≥1.75 in ≥80% of the samples.

Functional Testing (Fragment Purification and Restriction Digestion) - A linear pLA3 fragment is added to 8 Monarch[®] DNA Cleanup Columns, eluted and digested with DraIII-HF[®] resulting in ≥90% conversion to two-fragment pattern observed by gel electrophoresis in ≥75% of the samples.



Date 03 Dec 2017

Derek Robinson
Director of Quality Control

