

New England Biolabs Product Specification

Product Name: Monarch[®] RNA Cleanup Kit (10 µg)
Catalog #: T2030S/L
Kit Components: Monarch[®] RNA Cleanup Columns (10 µg) (T2037)
Monarch[®] Collection Tubes II (T2018)
Monarch[®] RNA Cleanup Binding Buffer (T2041)
Monarch[®] RNA Cleanup Wash Buffer (T2042)
Nuclease-free Water (B1500)

Shelf Life: 18 months
Storage Temp: 25°C
Specification Version: PS-T2030S/L v1.0
Effective Date: 25 May 2018

Assay Name/Specification (minimum release criteria)

Functional Testing (Enzyme Inactivation and Removal) - Eight Monarch[®] RNA Cleanup Columns (T2037, 10 µg) are tested with 1 µg of rRNA and 2 units of DNase I. A 50 µl reaction in DNase I Reaction Buffer with 1 µg of pBR322 and 50% of the column eluate incubated for 10 minutes at 37°C results in >99.95% removal of DNase I activity in ≥75% of the samples.

Functional Testing (RNA Recovery and Purity) - Twenty-four Monarch[®] RNA Cleanup Columns (T2037, 10 µg) are tested with 1 µg of rRNA resulting in ≥75% recovery in ≥80% of the samples. OD 260/280 and 260/230 are ≥1.75 in ≥80% of the samples.

* **Individual Product Component Note** - Standard Quality Control Tests are performed for each component included in Monarch[®] RNA Cleanup Kit (10 µg) and meet the designated specifications.



Date 25 May 2018

Derek Robinson
Director of Quality Control

