
New England Biolabs Product Specification

Product Name: Monarch[®] RNA Cleanup Columns (10 µg)
Catalog #: T2037L
Shelf Life: 18 months
Storage Temp: 25°C
Specification Version: PS-T2037L v1.0
Effective Date: 17 May 2018

Assay Name/Specification (minimum release criteria)

Functional Testing (RNA Recovery and Purity) - Twenty-four Monarch [®] RNA Cleanup Columns (T2037, 10 µg) are tested with 1 µg of rRNA resulting in ≥75% recovery in ≥80% of the samples. OD 260/280 and 260/230 are ≥1.75 in ≥80% of the samples.
--



Date 17 May 2018

Derek Robinson
Director of Quality Control

