240 County Road Ipswich, MA 01938-2723 Tel 978-927-5054 Fax 978-921-1350 www.neb.com info@neb.com

New England Biolabs Product Specification

Product Name: HiScribe® T7 High Yield RNA Synthesis Kit

Catalog #: E2040S

Kit Components: ATP (N0451)

GTP (N0452) UTP (N0453) CTP (N0454)

10X T7 Reaction Buffer (B2041) FLuc Control Template (N0426) T7 RNA Polymerase Mix (M0255)

0.1M DTT (B1222)

Shelf Life: 18 months
Storage Temp: -20°C

Specification Version: PS-E2040S v2.0 Effective Date: 15 Aug 2023

Assay Name/Specification (minimum release criteria)

Functional Testing (RNA Synthesis, Yield) - A 20 μ l reaction in 1X T7 Reaction Buffer in the presence of 40 mM NTPs containing 1 μ g FLuc Control Template with 2 μ l of the T7 RNA Polymerase Mix incubated for 2 hours at 37°C results in a minimum of 150 μ g of a 1.8 kb RNA as determined by OD 260 and agarose gel electrophoresis.

* Individual Product Component Note - Standard Quality Control Tests are performed for each component included in HiScribe® T7 High Yield RNA Synthesis Kit and meet the designated specifications.

One or more products referenced in this document may be covered by a 3rd-party trademark. Please visit www.neb.com/trademarks for additional information.

Lauren Brown Quality Approver







www.Brown

Date 15 Aug 2023