240 County Road Ipswich, MA 01938-2723 Tel 978-927-5054 Fax 978-921-1350 www.neb.com info@neb.com

New England Biolabs Product Specification

Product Name: HiScribe® T7 Quick High Yield RNA Synthesis Kit

Catalog #: E2050L

Kit Components: T7 RNA Polymerase Mix (M0255)

NTP Buffer Mix (N2052)
FLuc Control Template (N0426)
DNase I (RNase-free) (M0303)

LiCl Solution (B2051) 0.1M DTT (B1222)

Shelf Life: 24 months
Storage Temp: -20°C

Specification Version: PS-E2050L v1.0
Effective Date: 12 Feb 2024

Assay Name/Specification (minimum release criteria)

Functional Testing (RNA Synthesis, Yield) - A 20 μ l reaction in the presence of 40 mM NTPs containing 1 μ g FLuc Control Template with 2 μ l of the T7 RNA Polymerase Mix incubated for 2 hours at 37°C results in a minimum of 150 μ g of a 1.8 kb RNA as determined by OD 260 and agarose gel electrophoresis.

* Individual Product Component Note - Standard Quality Control Tests are performed for each component included in HiScribe® T7 Quick High Yield RNA Synthesis Kit and meet the designated specifications.

One or more products referenced in this document may be covered by a 3rd-party trademark. Please visit www.neb.com/trademarks for additional information.

Lauren Brown Quality Approver







aurenbrow

Date

12 Feb 2024