
New England Biolabs Product Specification

Product Name: Anti-CBD Monoclonal Antibody
Catalog #: E8034S
Recommended Dilution: 1:1,000
Shelf Life: 24 months
Storage Temp: -20°C
Storage Conditions: 150 mM NaCl, 10 mM HEPES, 50 % Glycerol, (pH 7.5 @ 25°C)
Specification Version: PS-E8034S v1.0
Effective Date: 14 Feb 2018

Assay Name/Specification (minimum release criteria)
Functional Testing (Western Blot, Antibody) - A 1:1,000 dilution of the antibody is tested for its ability to detect 50 ng of CBD fusion protein by Western Blot.



Date 14 Feb 2018

Derek Robinson
Director of Quality Control

