

## New England Biolabs Product Specification

*Product Name:* Monarch<sup>®</sup> RBC Lysis Buffer  
*Catalog #:* T3051L  
*Shelf Life:* 24 months  
*Storage Temp:* 4°C  
*Composition:* Proprietary  
*Specification Version:* PS-T3051L v1.0  
*Effective Date:* 16 Dec 2020

### Assay Name/Specification (minimum release criteria)

**Functional Testing (Genomic DNA Yield and Integrity Analysis)** - Genomic DNA is purified from 8 samples of whole pig blood using the High Molecular Weight DNA Extraction from Blood (Frozen Mammalian Blood) protocol. Each set of reagents is functionally tested to ensure successful isolation of gDNA from frozen pig blood samples. Yield of genomic DNA is compared across various metrics including a minimum gDNA yield of 6 µg per extraction and a recovery equivalence of ≥85% as compared to the control lot. OD 260/280 and 260/230 are ≥1.75 in ≥80% of the samples, and DIN values are ≥9 in ≥80% of the samples. In ≥ 80% of samples, ≥ 75% of the DNA is 50 kb in size or greater.

One or more products referenced in this document may be covered by a 3rd-party trademark.  
Please visit [www.neb.com/trademarks](http://www.neb.com/trademarks) for additional information.



Date 16 Dec 2020

Derek Robinson  
Director, Quality Control

