240 County Road Ipswich, MA 01938-2723 Tel 978-927-5054 Fax 978-921-1350 www.neb.com info@neb.com

New England Biolabs Product Specification

Product Name: pTXB1 Vector

Catalog #: N6707SConcentration: $200 \mu g/ml$ Shelf Life: 36 monthsStorage Temp: $-20^{\circ}C$

Storage Conditions: 10 mM Tris-HCl, 1 mM EDTA, (pH 8.0 @, 25°C)

Specification Version: PS-N6707S v1.0
Effective Date: 17 Jul 2018

Assay Name/Specification (minimum release criteria)

A260/A280 Assay - The ratio of UV absorption of pTXB1 Vector at 260 and 280 nm is between 1.8 and 2.0.

DNA Concentration (A260) - The concentration of pTXB1 Vector is between 200 and 220 μ g/ml as determined by UV absorption at 260 nm.

Restriction Digest (Linearization) - A 50 μ l reaction in CutSmart® Buffer containing 2 μ g of pTXB1 Vector and 10 units of SapI incubated for 1 hour at 37°C produces > 95% linearization resulting in a single band of approximately 6706 bp as determined by agarose gel electrophoresis.

Kuk Cotum

Date 17 Jul 2018

Derek Robinson

Director of Quality Control





