
New England Biolabs Product Specification

Product Name: Oligo d(T)₁₈ mRNA Primer
Catalog #: S1316S
Shelf Life: 36 months
Storage Temp: -20°C
Composition: Supplied as dried, triethylammonium salt
Specification Version: PS-S1316S v2.0
Effective Date: 22 Nov 2021

Assay Name/Specification (minimum release criteria)

Electrophoretic Pattern (Oligonucleotide) - The banding pattern of Oligo d(T) ₁₈ mRNA Primer on a 15% TBE-Urea gel is evaluated against a control lot for expected size, sharpness, and relative intensity.
--

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 22 Nov 2021

Derek Robinson
Director, Quality Control

